



Mirahold and Spandex Instructions For Use (IFU)

Limitation of reprocessing: max. 50 x

Place of handling: Clean surface with a disposable cloth or paper towel.

Storage and transport: Storage and Transport.

Preparations for cleaning: n/a

Cleaning: Only use cleansers suitable for synthetic material, carefully following the instructions of the respective manufacturer. After cleaning, remove cleanser thoroughly under running water. Do not exceed 200° F (93° C) for mechanical cleaning.

Disinfection: Only use disinfecting solution suitable for synthetic material, carefully following of the instructions the respective manufacturer. Remove disinfecting solution thoroughly under running water. Do not exceed 200° F (93° C) for mechanical disinfection.

Maintenance: No special requirements.

Packaging: Standardized packing material for sterilization can be used (recommendation).

Sterilization: Vapor sterilization at 250° F (121° C) for 15 minutes.

Control/Functional Check: Sight check on damages, wear, deformation.

Storage: No special requirements.

The a.m. instructions have been validated as **SUITABLE** for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine control are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.