## EN Uni-Syringe/Durasyringe Cleaning

Clean the external surfaces of the syringe using a solution of mild detergent and warm water. A soft bristled brush may be used to clean around the buttons, tip collar and syringe threads. Do not use powdered cleansers, scouring pads or abrasive scrubbers any of which can damage the surface finish of the syringe.

### Uni-Syringe/Durasyringe Disinfecting

The use of an EPA-registered hospital disinfectant is adequate. Follow disinfectant manufacturer's instructions for use when disinfecting syringe heads.

### Syringe Tip Cleaning & Sterilization

Clean tips under running water for 30 seconds with a soft brush and place in an ultrasonic bath with an enzymatic cleaner to remove superficial debris prior to sterilization; any debris may reduce the effectiveness of the sterilization.

Sterilization process is by steam autoclave. Parameters for wrapped are 132°C (270°F) for 15 minutes at temperature with a 30-minute drying time when using a gravity displacement autoclave. Parameters for unwrapped are 132°C (270°F) for 3 minutes at temperature when using a gravity displacement autoclave.

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