

Now that you can sterilize your suction handpieces after each patient... **why take any chances?**



Long HVE Medium HVE Short HVE SS S.E. S.E.



Item # Description

Bull Frog "Expert" Series H.V.E.

B11190	Bull Frog "Expert" Series - Long HVE - Includes 1/2" (12.7 mm) Swivel Tubing Adapter (1) (100% Stainless Steel)
B11200	Bull Frog "Expert" Series - Medium HVE - Includes 1/2" (12.7 mm) Swivel Tubing Adapter (1) (100% Stainless Steel)
B11195	Bull Frog "Expert" Series - Short HVE - Includes 1/2" (12.7 mm) Swivel Tubing Adapter (1) (100% Stainless Steel)

Bull Frog "Expert" Series S.E.

B11130	Bull Frog "Expert" Series S.S. (100% Stainless Steel) - Includes Swivel Tubing Adapter & Gray Silicone Tip
B22270	Bull Frog "Expert" Series S.E. (100% Stainless Steel) - Includes Swivel Tubing Adapter

Introducing the BullFrog "Expert" Series Suction Handpieces

Unlike less durable aluminum versions, the BullFrog "Expert" Series is a 100% stainless steel. This means they are completely autoclavable and withstand sterilization after each and every patient.

- ✓ Available in High Volume and Saliva Ejectors
- ✓ Safe to use with sodium hypochlorite
- ✓ Meets the U.S. government's agency mandate for autoclaving after every patient
- ✓ 5 Year Warranty

In addition, the BullFrog "Expert" Series features a sleek ergonomic body for better grip and balance as well as a newly designed metal lever for smooth opening and closing action.

**5 Year
Warranty**

Did You Know?

The CDC recommends sterilizing suction handpieces after every patient and it's a requirement for VA and Military dental facilities.



(800) 328-2335 | www.hagerworldwide.com

Now that you have selected the Bull Frog "Expert" Series... Compliment with autoclavable accessories



HAGER "DENTAL INSTRUMENT" CASSETTES

Item #	Description	Qty.
B1501	HAGER "Dental Instrument" Cassettes - for Five (5) Instruments Dimensions: 8" x 3" x 1"/20.3 cm x 7.6 cm x 2.5 cm (100% Stainless Steel)	1
B1502	HAGER "Dental Instrument" Cassettes - for Eight (8) Instruments Dimensions: 8" x 5-1/2" x 1"/20.3 cm x 14 cm x 2.5 cm (100% Stainless Steel)	1
B1503	HAGER "Dental Instrument" Cassettes - for Oral Surgery Dimensions: 11" x 8" x 1"/27.9 cm x 20.3 cm x 2.5 cm (100% Stainless Steel)	1


VACUUM TIPS SS

Item #	Description	Qty.
B10100	Stainless Tip - 6" (15 cm) - Stainless Steel (Autoclavable)	1
B10102	Stainless Tip - 6" (15 cm) - Vented - Stainless Steel (Autoclavable)	1
B10101	Stainless Tip - 4-1/2" (11 cm) - Stainless Steel (Autoclavable)	1
B10103	Stainless Tip - 4-1/2" (11 cm) - Vented - Stainless Steel (Autoclavable)	1



PELOTTE PLUS EVACUATION TIPS

Item #	Description	Qty.
254024	Green - Length: 4" (10 cm) by 7/16" diameter (11 mm) (Autoclavable)	10/Box



SVEDOPTER OPTICAL SALIVA EJECTOR

Item #	Description	Qty.
254140	Saliva Ejector Optical w/2 Mirrors (Autoclavable)	1



SURGICAL TIP SURGICAL ASPIRATION TIP SS

Item #	Description	Qty.
B10104	Stainless Steel - Diameter: 3/16" (4.5 mm) Length 6" (15.2cm) (Autoclavable)	1

