


# A New Angle on Reducing Hand and Wrist Fatigue.



Ergonomic ingenuity enables the Ergo-Vac HVE handpiece to relieve awkward wrist and hand positioning while suctioning. With the tip at a 110° angle to the handpiece, you enjoy a more neutral hand/wrist posture and reduced hand and arm fatigue. Plus, Ergo-Vac provides the dentist with improved access to the work area. Ergo-Vac also features:

- One-handed operation thanks to convenient finger-tip lever
- 360° rotating coupling prevents hose from turning or kinking
- Quick disconnect for easy cleaning and autoclaving
- Swivel or stationary head

Ask your dealer for Ergo-Vac, or call 800.328.2335 for the dealer nearest you or additional information.



360° swivel head.

*Performance You Can Depend On™*



## Ergo•Vac

Tel: (800) 328-2335 [hagerworldwide.com](http://hagerworldwide.com)

Ergo-Vac is a registered trademark of Hager Worldwide, Inc.

# Ergo-Vac Can Minimize Repetitive Strain Injuries (RSI)

Repetitive Strain Injuries (RSIs) are injuries to muscles, tendons and nerves that become irritated or inflamed due to repetitive motions, excessive force, and extreme motions. And RSIs are on the rise in dentistry. You can minimize possible RSIs with the Ergo-Vac high volume evacuation handpiece (HVE). By positioning the wrist at 110° to the patient, versus the traditional 0°, Ergo-Vac helps you assume a healthier, more ergonomic posture. At the end of the day, you will have less fatigue, you will experience less pain, and you will be able to provide the quality of service that you and your patient demand.



## Ergo-Vac HVE Handpiece

- w/ Swivel Adapter (1) #B11171
- w/o Swivel Adapter (1) #B11181

## Ergo-Vac Combo Kit

- Contains: Ergo HVE Handpieces (3)
- 1/2" (1.2cm) Swivel Tubing Adapter (1) #B11180

## Ergo-Vac Stationary Head

- w/ Swivel Adapter (1) #B11171S
- w/o Swivel Adapter (1) #B11181S

## Ergo-Vac Stationary Combo Kit

- Contains: Ergo HVE Handpieces (3)
- 1/2" (1.2cm) Swivel Tubing Adapter (1) #B11180S



Optional Swivel Adapter



Stationary Combo Kit

# More Performance and Productivity Products from Hager Worldwide



## Polymatic

Keeps tongue and cheek away, while removing spray mist and secretions with ease!

- Dual-ended aspiration at same time
- Simultaneous cheek and tongue retraction
- Provides clear, free work area
- No suck-down of soft tissue
- Autoclavable plastic 340°F (171°C)

Can be used for any type of procedure with a high/low speed handpiece, air/water syringe or scaling.

**Polymatic**  
#250435 Set of 2 (left and right)



## Mirror Suction

Mouth Mirror with Mouth Suction - perfect for pit and fissure sealants and pedodontic treatment!

- Acts as a tongue guard
- Provides excellent aspiration
- Front surface mirror
- No assistant required

### "Mini" Kit

#254068 2 Handles w/ Mirrors,  
1 Hose, 1 11mm Adapter,  
1 Cleaning Brush

### Mirror Suction Introductory Kit

#254050 5 Handles w/ Mirrors,  
2 Hoses, 2 Adapters



## Vacuum Tips

Stainless steel high volume evacuation tips.

- Vented and non-vented (no suck-down of soft tissue)
- Seamless for easy cleaning
- Autoclavable

### Vacuum Tip Stainless Steel

#B10101 4 1/2" non-vented  
#B10103 4 1/2" vented  
#B10100 6" non-vented  
#B10102 6" vented