


Your Implant Impression Doesn't Have To Be Complicated.



Simplify it with Miratray Implant Advanced

- ✓ Clear, pierceable surface eliminates the need for cutting holes in a stock tray - saving you time and hassle
- ✓ Compatible with all implant systems and impression materials, so there is no need to purchase additional supplies



(800) 328-2335

www.hagerworldwide.com

What your fellow dentists are saying...

I have always preferred the open-tray technique for implant abutment transfer pick-up impressions, both single and multiple units. What I don't like, however, is having to custom-make a tray or cut open a window in a stock tray. Plus, there's always the concern that the impression material could come through and obliterate the screw pin holding the transfer abutment.

Hager Worldwide's Miratray Implant Advanced system eliminates all of those issues and makes it easy for you to take precise pick-up impressions.

Howard S. Glazer, DDS, FAGD

The Hager Miratray Implant Advanced tray simplifies open tray implant impressions by eliminating the issues with modifying stock trays. This assures that once filled with impression material and inserted, you can get the pins to come through the hole in the tray.

Gregori M. Kurtzman, DDS, MAGD, DICOI

Take our Free 1 hour CE course

*Dental Implant Impressions:
Steps to Easier and More
Predictable Results*

www.hagerworldwide.com/CEcourses



Item Number	Description	Quantity
# 101250	Introductory Kit	1 of each size
# 101251	Upper S1 (small)	Set of 6
# 101252	Upper S2 (medium)	Set of 6
# 101253	Upper S3 (large)	Set of 6
# 101254	Lower I1 (small)	Set of 6
# 101255	Lower I2 (medium)	Set of 6
# 101256	Lower I3 (large)	Set of 6



*Sold in packages of 6
Order through your preferred dealer today!*



(800) 328-2335

www.hagerworldwide.com