

HAGER WORLDWIDE, INC.

SAFETY DATA SHEET

SECTION 1-PRODUCT IDENTIFICATION

PRODUCT NAME MIRA WHITE OXYGEN
MANUFACTURER HAGER WORLDWIDE, INC.
441 19TH STREET SE
HICKORY, NC 28602
TEL: 800-328-2335
FAX: 800-573-9392

EMERGENCY NUMBER 800-328-2335

SECTION 2-HAZARDS IDENTIFICATION

Emergency Overview

Clear Gel
Slight eye irritation potential
Not an environmental hazard

POTENTIAL HEALTH EFFECTS

PRIMARY ENTRY ROUTES: ORAL, EYE

ACUTE EFFECTS:

INHALATION: NONE KNOWN

EYE: MAY CAUSE SLIGHT IRRITATION, BUT WILL NOT CAUSE INJURY

SKIN: NOT A SKIN IRRITANT, PROLONGED CONTACT MAY CAUSE
TRANSIENT IRRITATION.

INGESTION: NOT CONSIDERED TOXIC BY INGESTION. SWALLOWING SMALL
AMOUNTS (1 OZ. OR MORE) MAY CAUSE NAUSEA AND VOMITING
WITHIN MINUTES AND DIARRHEA AND MILD GI IRRITATION MAY
OCCUR UP TO 24 HOURS AFTER INGESTION.

MEDICAL CONDITIONS AGGRAVATED BY LONG-TERM EXPOSURE: NONE KNOWN

CHRONIC EFFECTS/CARCINOGENICITY: NOT APPLICABLE

SECTION 3-COMPOSITION/INFORMATION ON INGREDIENTS

INGREDIENT NAME	CAS NUMBER
WATER	7732-18-5
GLYCERIN	56-81-5
CARBAMID PEROXIDE	124-43-6
POLYETHYLENGLYCOL	25322-68-3

SECTION 4-FIRST AID MEASURES

EYE CONTACT:	RINSE THOROUGHLY WITH WATER, IF IRRITATION PERSISTS, SEEK MEDICAL ATTENTION.
SKIN CONTACT:	RINSE THOROUGHLY WITH WATER.
INGESTION:	FOR INGESTION OF SMALL AMOUNTS (1 OZ OR MORE), ADMINISTER MILK (WATER OR OTHER LIQUID IF MILK IS UNAVAILABLE) TO DILUTE STOMACH CONTENTS. FOR INGESTION OF LARGER AMOUNTS (>3 OZ. IN CHILDREN, > 8 OZ. IN ADULTS), ADMINISTER MILK (WATER OR OTHER LIQUID IF MILK IS UNAVAILABLE) TO DILUTE STOMACH CONTENTS AND CONTACT A PHYSICIAN FOR FURTHER TREATMENT.
INHALATION:	NOT APPLICABLE.

SECTION 5-FIRE-FIGHTING MEASURES

FLASH POINT:	NONE
FLASH POINT METHOD:	NOT APPLICABLE
LEL:	NOT APPLICABLE
UEL:	NOT APPLICABLE
FLAMMABILITY CLASSIFICATION:	NON-FLAMMABLE
EXTINGUISHING MEDIA:	WATER, DRY CHEMICAL, ALCOHOL FOAM, CO2
UNUSUAL FIRE OR EXPLOSION HAZARD:	NONE KNOWN
HAZARDOUS COMBUSTION PRODUCTS:	NONE KNOWN

FIRE-FIGHTING INSTRUCTIONS:

DO NOT RELEASE RUN-OFF FROM FIRE CONTROL METHODS TO SEWERS OR WATERWAYS.

FIRE-FIGHTING EQUIPMENT:

BECAUSE FIRE MAY PRODUCE TOXIC THERMAL DECOMPOSITION PRODUCTS, WEAR A SELF-CONTAINED BREATHING APPARATUS WITH A FULL FACEPIECE OPERATED IN PRESSURE-DEMAND OR POSITIVE-PRESSURE MODE.

SECTION 6-ACCIDENTAL RELEASE MEASURES

SCOOP UP SOLIDS INTO A WASTE CONTAINER. ABSORB LIQUIDS WITH APPROPRIATE INERT ABSORBENT. WASH AWAY SMALL QUANTITIES OF RESIDUAL MATERIAL WITH WATER. PRODUCT SPILLS MAY PRESENT A SLIPPING HAZARD. USE CAUTION DURING CLEAN-UP AND FLUSH AREA THOROUGHLY WITH WATER.

SECTION 7-HANDLING AND STORAGE

FOR BULK CONTAINERS AND CASE QUANTITIES:

KEEP CONTAINERS CLOSED AND STORE PRODUCT AT AMBIENT TEMPERATURES AWAY FROM DIRECT SUNLIGHT. AVOID FREEZING OR EXCESSIVE HEAT. AVOID CONTACTING PRODUCT WITH EYES.

SECTION 8-EXPOSURE CONTROLS/PERSONAL PROTECTION

RESPIRATORY PROTECTION:

NONE REQUIRED

PROTECTIVE GLOVES:

NOT NORMALLY REQUIRED. WEAR GENERAL-PURPOSE GLOVES IF HANDLING BULK PRODUCT.

EYE PROTECTION:

NOT NORMALLY REQUIRED. WEAR SAFETY GLASSES OR GOGGLES IF HANDLING BULK PRODUCT.

OTHER PROTECTIVE CLOTHING/EQUIPMENT:

CLEAN BODY COVERING CLOTHING IS SUFFICIENT

SECTION 9-PHYSICAL AND CHEMICAL PROPERTIES

APPEARANCE: TRANSLUCENT GEL
ODOR: NONE
ODOR THRESHOLD: NOT KNOWN
VAPOR PRESSURE: NOT DETERMINED
VAPOR DENSITY (Air=1): NOT DETERMINED
SPECIFIC GRAVITY (H₂O=1, AT 4 DEGREES CELSIUS): 1.25 g/cm³
Ph: ca. 6-7
WATER SOLUBILITY: MODERATELY SOLUBLE
BOILING POINT: 100 DEGREES CELSIUS
FREEZING/MELTING POINT: 5 DEGREES CELSIUS
EVAPORATION POINT: NOT KNOWN

SECTION 10-STABILITY AND REACTIVITY

STABILITY: PRODUCT IS STABLE AT ROOM TEMPERATURE IN CLOSED CONTAINERS UNDER NORMAL STORAGE AND HANDLING CONDITIONS.
POLYMERIZATION: HAZARDOUS POLYMERIZATION WILL NOT OCCUR.
CHEMICAL INCOMPATIBILITIES: NONE KNOWN
CONDITIONS TO AVOID: NONE KNOWN
HAZARDOUS DECOMPOSITION PRODUCTS: NONE KNOWN

SECTION 11-TOXICOLOGICAL INFORMATION

TOXICITY DATA:

EYE EFFECTS: SLIGHT TRANSIENT IRRITATION POSSIBLE
SKIN EFFECTS: NOT A PRIMARY IRRITANT
ACUTE INHALATION EFFECTS: NOT APPLICABLE
ACUTE ORAL EFFECTS: CONSIDERED NON-TOXIC
CHRONIC EFFECTS: NONE KNOWN
CARCINOGENICITY: NO
MUTAGENICITY: NO
TERATOGENICITY: NO

SECTION 12-ECOLOGICAL INFORMATION

THIS PRODUCT IS NOT CLASSIFIED AS ENVIRONMENTALLY HAZARDOUS.

SECTION 13-DISPOSAL CONSIDERATIONS

ANY DISPOSAL PRACTICE MUST BE IN COMPLIANCE WITH LOCAL, STATE, AND FEDERAL REGULATIONS. LOCAL AND STATE LAWS MAY BE DIFFERENT THAN FEDERAL ONES. BE SURE TO CONSULT WITH LOCAL AND STATE AGENCIES FOR SPECIFIC RULES. LANDFILL THIS PRODUCT IN ACCORDANCE WITH ALL APPLICABLE REGULATIONS. DILUTE AND DRAIN SMALL AMOUNTS TO SEWER.

SECTION 14-TRANSPORTATION INFORMATION

D.O.T. SHIPPING NAME:	NON-HAZARDOUS
TECHNICAL SHIPPING NAME:	TOOTH WHITENING GEL
D.O.T. HAZARD CLASS:	NONE
U.N./N.A. NUMBER:	NONE
HAZARDOUS SUBSTANCE/RQ:	NONE
D.O.T. LABEL:	NONE
D.O.T. PLACARD:	NONE

SECTION 15-REGULATORY INFORMATION

MIRA WHITE OXYGEN IS A COSMETIC PRODUCT THAT IS REGULATED BY THE UNITED STATES FOOD AND DRUG ADMINISTRATION.

SECTION 16-OTHER INFORMATION

DISCLAIMER: THIS SAFETY DATA SHEET IS OFFERED SOLELY FOR YOUR INFORMATION AND CONSIDERATION. THE INFORMATION PRESENTED HERE IS ACCURATE AND RELIABLE TO THE BEST OF OUR KNOWLEDGE. HOWEVER, NO WARRANTIES (EITHER EXPRESSED OR IMPLIED) ARE PROVIDED AND WE ASSUME NO RESPONSIBILITY FOR THE ACCURACY OR COMPLETENESS OF DATA CONTAINED HEREIN. PERSONS RECEIVING THIS INFORMATION ARE URGED TO MAKE THEIR OWN DETERMINATION AS TO THE INFORMATION SUITABILITY FOR THEIR PARTICULAR APPLICATION.

ISSUING DEPARTMENT:
TELEPHONE:

HAGER WORLDWIDE, INC.
800-328-2335

Revision Date: 11/20/2015

Printed Date: 11/20/2015