

Hygienists, meet your new assistant...

The Mirror Suction

Imagine having your mouth mirror, aspiration, and retraction all in one instrument.

Water flow can be better controlled during ultra-sonic scaling. Sealants can be applied without assistance.

That's why the Mirror Suction provides strong suction power by connecting to your HVE valve.

- **Suit Your Preference** with double and single-sided mirrors
- **Easy-to-Clean** because all components are autoclavable
- **Reduces Fogging** through airflow in the mirror head
- **For best results** use in tandem with the saliva ejector

Order from your dealer using the below part numbers.



Double-Sided (#5 mirror head)

Mirror Suction - Double-Sided

254069	Mirror Suction Double-Sided "Introductory" Kit: (5) Handles with #5 double-sided mirrors, (2) Hoses, (2) 11 mm adapters and (1) Cleaning brush
254070	Mirror Suction Double-Sided - "Mini" Kit - (2) Handles with (2) #5 double-sided mirrors, (1) Hose, (1) 11 mm adapter & (1) Cleaning brush all packaged in a sterilizable cassette
254077	Mirror Suction Double-Sided #5 Mirror Head - Pack of 6 heads



Single-Sided (#4 mirror head)

Mirror Suction - Single-Sided

254050	Mirror Suction "Introductory" Kit: (5) Handles with #4 mirrors, (2) Hoses, (2) 11 mm adapters and (1) Cleaning brush
254068	Mirror Suction "Mini" Kit: (2) Handles with #4 mirrors, (1) Hose, (1) 11 mm adapter and (1) Cleaning brush
254058	Mirror Suction "Single-Sided" #4 Mirror Head - Pack of 6 heads



Mirror Suction Introductory Kit



Mirror Suction Mini Kit



(800) 328-2335

www.hagerworldwide.com